





Holder Type	Description	Applications	Maximum Size mm (in) WxHxD	Exposure Size mm (in) WxH	Instrument	Rack Capacity
 SL-3T Single exposure window with spring clip back Part Number 19163900	General purpose holder for exposure of thin specimens	Textiles, plastic film, automotive interior	67 x 145 x 3 (2.6 x 5.7 x 0.1)	50 x 121 (2.0 x 4.7)	Ci3000+ Ci4000 Ci5000	20 68 111
 SL-3T with Glass Single exposure window with glass and adjustable back Part Number 07303900	Exposure through secondary filter to closely simulate end-use environment	Textiles, paper, plastic film, carpet, automotive interior	67 x 145 x 15 (2.6 x 5.7 x 0.6)	50 x 121 (2.0 x 4.7)	Ci3000+ Ci4000 Ci5000	20 68 111
 CD-3T Three exposure windows with spring clip back Part Number 20215700	General purpose holder with three exposure windows for thin specimens	Textiles, paper, plastic film, automotive interior	67 x 145 x 3 (2.6 x 5.7 x 0.1)	3 windows: 38 x 50 (1.5 x 2.0)	Ci3000+ Ci4000 Ci5000	20 68 111
 CD-3T with Glass Three exposure windows with glass, spring clip back Part Number 07303800	Exposure through secondary filter to closely simulate end-use environment	Textiles, paper, plastic film, wood, automotive interior	67 x 145 x 15 (2.6 x 5.7 x 0.6)	3 windows: 38 x 50 (1.5 x 2.0)	Ci3000+ Ci4000 Ci5000	20 68 111
 RD-3T Single window or 3-window exposure with "bulldog" clip holder Part Number 20017900	Versatile holder to accommodate a range of specimen widths and thicknesses	Coatings on various substrates, plastics, textiles, glass	77 x 152 x 10 (3.0 x 6.0 x 0.4)	57 x 134 (2.3 x 5.3)	Ci3000+ Ci4000 Ci5000	20 68 111
 TEX-3T with Mask Single exposure window with mask, adjustable Part Number 19186700	Holder with masked exposure and adjustable back to allow for thick sample	Textiles, foam, foam-backed materials	45 x 134 x 12 (1.7 x 5.3 x 0.5)	19 x 119 (0.7 x 4.7)	Ci3000+ Ci4000 Ci5000	29 104 170
 Panel Holder Part Number 19188501	Holds 76 x 152 (3 x 6) panels for exposure tests of coatings	Coatings, rigid plastic, wood	76 x 152 x 9 (3.0 x 6.0 x 0.37)	76 x 146 (3.0 x 5.7)	Ci3000+ Ci4000 Ci5000	16 56 87
 Drop-in Specimen Bar Holder Part Number 19184600	Holds vertical specimen bars	Plastics	77 x 144 x 3 (3.0 x 5.6 x 0.1)	76 x 125 (3.0 x 4.9)	Ci3000+ Ci4000 Ci5000	15 50 84

This chart is a representative sample of specimen holders available for a variety of Atlas® instruments. For specific information about specimen holders that best meet your needs, please contact your local Atlas representative.



# Specimen Holders

## Accelerated Weathering Products

Rack Capacity	Instrument	Exposure Size mm (in) WxH	Maximum Size mm (in) WxHxD	Applications	Description	Holder Type
11	Xenotest® Alpha+ Xenotest 150 S+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Textiles, plastics, coatings, paper	Holds specimen up to 3 mm thick	Regular Specimen Holder Part Number 56075142
11	Xenotest Alpha+ Xenotest 150 S+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Carpets, foam, foam-backed materials, thick panels	Holds specimen up to 15 mm thick	Specimen Holder for Thick Specimens Part Number 56077987
only 1 needed (replaces one holder)	Xenotest Alpha+ Xenotest 150 S+	—	135 x 45 (5.3 x 1.8)	Blue Scale fabric	Holder for Blue Scale fabric during weathering tests	Specimen Holder for Blue Scale Part Number 56050873
16	Xenotest Beta+	295 x 70 (11.6 x 2.8)	310 x 80 (12.2 x 3.2)	Textiles, plastics, coatings, paper, window profiles	Holds specimen up to 10 mm thick	Standard Specimen Holder (including back plate) Part Number 56076543
16	Xenotest Beta+	90 x 58 (3.5 x 2.3)	100 x 68 (3.9 x 2.7)	Textiles, plastics, coatings, paper	Holds 3 specimens up to 10 mm thick, fixed by a back plate and a cover mask	Special Specimen Holder, 3 Segments Part Number 56076699
28 19	Xenotest Beta+ Xenotest 220/220+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Carpet, foam- backed textiles	Holds 2 specimens up to 10 mm thick, fixed by 2 cover masks	Special Specimen Holder, 2 Segments Part Number 56076700
28 19	Xenotest Beta+ Xenotest 220/220+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Textiles, paper, leather	Holds 2 samples up to 3 mm thick, fixed by 2 cover masks	Specimen Holder Standard, 2 Segments Quick-insertion Part Number 56079776
28 19	Xenotest Beta+ Xenotest 220/220+	270 x 45 (10.6 x 1.8)	320 x 48 (12.6 x 1.9)	Technical textiles	Holds 1 specimen up to 3 mm thick	Specimen Holder Special TS, 1 Segment Quick-insertion, Part Number 56079797
only 1 needed (replaces one holder)	Xenotest 220/220+	—	135 x 45 (5.3 x 1.8)	Blue Scale fabric	Holds Blue Scale fabric for weathering tests	Specimen Holder for Blue Scale (Special holder kit, 2 Segments) Part Number 56078036
15	SUNTEST® XLS+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Plastics, coatings, paper, leather	General purpose holder for various materials	Set of 15 Specimen Trays Part Number 56079636
34	SUNTEST XXL/XXL+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Plastics, coatings, paper	General purpose holder for various materials	Set of Five Specimen Trays for 7 Specimens Part Number 56079259
8	SUNTEST XXL/XXL+	295 x 70 (11.6 x 2.8)	310 x 80 (12.2 x 3.2)	Plastics, coatings, paper	General purpose holder for various materials	Set of Two Specimen Trays for 4 Specimens Part Number 56079261



This chart is a representative sample of specimen holders available for a variety of Atlas® instruments. For specific information about specimen holders that best meet your needs, please contact your local Atlas representative.